



CD70 Antibody, FITC conjugated

Product Code	CSB-PA08009C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P32970
Immunogen	Recombinant Human CD70 antigen protein (62-185AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	CD70 antigen (CD27 ligand) (CD27-L) (Tumor necrosis factor ligand superfamily member 7) (CD antigen CD70), CD70, CD27L CD27LG TNFSF7
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CD70
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.