



TUBA4A Antibody, HRP conjugated

Product Code	CSB-PA09499B0Rb
Abbreviation	Tubulin alpha-4A chain
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P68366
Immunogen	Recombinant Human Tubulin alpha-4A chain protein (20-183AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Relevance	Tubulin is the major constituent of microtubules. It binds two moles of GTP, one at an exchangeable site on the beta chain and one at a non-exchangeable site on the alpha chain.
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tubulin alpha-4A chain (Alpha-tubulin 1) (Testis-specific alpha-tubulin) (Tubulin H2-alpha) (Tubulin alpha-1 chain), TUBA4A, TUBA1
Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	TUBA4A