



TUBA4A Antibody, Biotin conjugated

Product Code	CSB-PA09499D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P68366
Immunogen	Recombinant Human Tubulin alpha-4A chain protein (20-183AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tubulin alpha-4A chain (Alpha-tubulin 1) (Testis-specific alpha-tubulin) (Tubulin H2-alpha) (Tubulin alpha-1 chain), TUBA4A, TUBA1
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	TUBA4A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.