





NRG1 Antibody

Product Code	CSB-PA09589A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q02297
Immunogen	Recombinant Human Pro-neuregulin-1, membrane-bound isoform protein (75-176AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pro-neuregulin-1, membrane-bound isoform (Pro-NRG1) [Cleaved into: Neuregulin-1 (Acetylcholine receptor-inducing activity) (ARIA) (Breast cancer cell differentiation factor p45) (Glial growth factor) (Heregulin) (HRG) (Neu differentiation factor) (Sensory and motor neuron-derived factor)], NRG1, GGF HGL HRGA NDF SMDF
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NRG1

Image



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA09589A0Rb at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA09589A0Rb at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.