



NRG1 Antibody, HRP conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA09589B0Rb |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q02297 |
| Immunogen | Recombinant Human Pro-neuregulin-1, membrane-bound isoform protein (75-176AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Pro-neuregulin-1, membrane-bound isoform (Pro-NRG1) [Cleaved into: Neuregulin-1 (Acetylcholine receptor-inducing activity) (ARIA) (Breast cancer cell differentiation factor p45) (Glial growth factor) (Heregulin) (HRG) (Neu differentiation factor) (Sensory and motor neuron-derived factor)], NRG1, GGF HGL HRGA NDF SMDF |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Neuroscience |
| Target Names | NRG1 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |