



MAPT Antibody, Biotin conjugated

Product Code	CSB-PA11169D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10636
Immunogen	Recombinant Human Microtubule-associated protein tau protein (33-166AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Microtubule-associated protein tau (Neurofibrillary tangle protein) (Paired helical filament-tau) (PHF-tau), MAPT, MAPTL MTBT1 TAU
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	MAPT
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.