



AREG Antibody, FITC conjugated

Product Code	CSB-PA11339C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15514
Immunogen	Recombinant Human Amphiregulin protein (107-148AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Amphiregulin (AR) (Colorectum cell-derived growth factor) (CRDGF), AREG, AREGB SDGF
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	AREG
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.