



# AREG Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA11339C0Rb
<b>Abbreviation</b>	Amphiregulin
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P15514
<b>Immunogen</b>	Recombinant Human Amphiregulin protein (107-148AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Relevance</b>	Ligand of the EGF receptor/EGFR. Autocrine growth factor as well as a mitogen for a broad range of target cells including astrocytes, Schwann cells and fibroblasts.
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Amphiregulin (AR) (Colorectum cell-derived growth factor) (CRDGF), AREG, AREGB SDGF
<b>Species</b>	Human
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	AREG