




REN Antibody

Product Code	CSB-PA12069A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00797
Immunogen	Recombinant Human Renin protein (229-398AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IF; Recommended dilution: WB:1:2000-1:5000, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Renin (EC 3.4.23.15) (Angiotensinogenase), REN
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	REN

Image



Western Blot

Positive WB detected in: Mouse heart tissue, Mouse spleen tissue, Mouse lung tissue, Mouse stomach tissue


All lanes: REN antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 46, 45 kDa

Observed band size: 46, 29 kDa



Immunofluorescent analysis of HepG2 cells using CSB-PA12069A0Rb at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

