




REN Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA12069A0Rb |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P00797 |
| Immunogen | Recombinant Human Renin protein (229-398AA) |
| Raised In | Rabbit |
| Species Reactivity | Human, Mouse |
| Tested Applications | ELISA, WB, IF; Recommended dilution: WB:1:2000-1:5000, IF:1:50-1:200 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Renin (EC 3.4.23.15) (Angiotensinogenase), REN |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cardiovascular |
| Target Names | REN |

Image



Western Blot

Positive WB detected in: Mouse heart tissue, Mouse spleen tissue, Mouse lung tissue, Mouse stomach tissue


All lanes: REN antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 46, 45 kDa

Observed band size: 46, 29 kDa



Immunofluorescent analysis of HepG2 cells using CSB-PA12069A0Rb at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.