



SLC2A2 Antibody, FITC conjugated

Product Code	CSB-PA13329C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11168
Immunogen	Recombinant Human Solute carrier family 2, facilitated glucose transporter member 2 protein (67-174AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 2, facilitated glucose transporter member 2 (Glucose transporter type 2, liver) (GLUT-2), SLC2A2, GLUT2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC2A2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.