



AQP3 Antibody, Biotin conjugated

Product Code	CSB-PA13909D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92482
Immunogen	Recombinant Human Aquaporin-3 protein (210-244AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Aquaporin-3 (AQP-3) (Aquaglyceroporin-3), AQP3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	AQP3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.