



PRKCA Antibody, Biotin conjugated

Product Code	CSB-PA15859D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17252
Immunogen	Recombinant Human Protein kinase C alpha type protein (495-659AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein kinase C alpha type (PKC-A) (PKC-alpha) (EC 2.7.11.13), PRKCA, PKCA PRKACA
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PRKCA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.