



MAPK8 Antibody, HRP conjugated

Product Code	CSB-PA16759B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P45983
Immunogen	Recombinant Human Mitogen-activated protein kinase 8 protein (5-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitogen-activated protein kinase 8 (MAP kinase 8) (MAPK 8) (EC 2.7.11.24) (JNK-46) (Stress-activated protein kinase 1c) (SAPK1c) (Stress-activated protein kinase JNK1) (c-Jun N-terminal kinase 1), MAPK8, JNK1 PRKM8 SAPK1 SAPK1C
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MAPK8