



MAPK8 Antibody, HRP conjugated

Product Code	CSB-PA16759B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P45983
Immunogen	Recombinant Human Mitogen-activated protein kinase 8 protein (5-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitogen-activated protein kinase 8 (MAP kinase 8) (MAPK 8) (EC 2.7.11.24) (JNK-46) (Stress-activated protein kinase 1c) (SAPK1c) (Stress-activated protein kinase JNK1) (c-Jun N-terminal kinase 1), MAPK8, JNK1 PRKM8 SAPK1 SAPK1C
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MAPK8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.