



ATP1A1 Antibody, FITC conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA18129C0Rb |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P05023 |
| Immunogen | Recombinant Human Sodium/potassium-transporting ATPase subunit alpha-1 protein (225-407AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Sodium/potassium-transporting ATPase subunit alpha-1 (Na(+)/K(+)) ATPase alpha-1 subunit (EC 3.6.3.9) (Sodium pump subunit alpha-1), ATP1A1 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Signal Transduction |
| Target Names | ATP1A1 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |