



SRGAP2B Antibody, FITC conjugated

Product Code	CSB-PA20639C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0DMP2
Immunogen	Recombinant Human SLIT-ROBO Rho GTPase-activating protein 2B protein (79-150AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SLIT-ROBO Rho GTPase-activating protein 2B (SLIT-ROBO Rho GTPase activating protein 2 pseudogene 2), SRGAP2B, SRGAP2P2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SRGAP2B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.