




MEIKIN Antibody

Product Code	CSB-PA20649A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A0A087WXM9
Immunogen	Recombinant Human Meiosis-specific kinetochore protein (18-91AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Meiosis-specific kinetochore protein, MEIKIN
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MEIKIN

Image



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA20649A0Rb at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.