



PET100 Antibody, Biotin conjugated

Product Code	CSB-PA319771LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0DJ07
Immunogen	Recombinant Human Protein PET100 homolog, mitochondrial protein (26-73AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein PET100 homolog, mitochondrial, PET100, C19orf79
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PET100
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.