



Gapdh1 Antibody, Biotin conjugated

Product Code	CSB-PA362145LD01DLU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07486
Immunogen	Recombinant Drosophila melanogaster Glyceraldehyde-3-phosphate dehydrogenase 1 protein (1-332AA)
Raised In	Rabbit
Species Reactivity	Drosophila melanogaster
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glyceraldehyde-3-phosphate dehydrogenase 1 (EC 1.2.1.12) (Glyceraldehyde-3-phosphate dehydrogenase I) (GAPDH I), Gapdh1, Gadph-1
Immunogen Species	Drosophila melanogaster
Research Area	Others
Target Names	Gapdh1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.