




MEI4 Antibody

Product Code	CSB-PA427046LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A8MW99
Immunogen	Recombinant Human Meiosis-specific protein MEI4-like protein (22-137AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Meiosis-specific protein MEI4, MEI4
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MEI4

Image



Western Blot

Positive WB detected in: Mouse heart tissue

All lanes: MEI4 antibody at 2.7µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 45 kDa

Observed band size: 45 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.