



UNC79 Antibody, Biotin conjugated

Product Code	CSB-PA603875LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P2D8
Immunogen	Recombinant Human Protein unc-79 homolog protein (1562-1678AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein unc-79 homolog, UNC79, KIAA1409
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	UNC79