



KLF5 Antibody, HRP conjugated

Product Code	CSB-PA613429LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13887
Immunogen	Recombinant Human Krueppel-like factor 5 protein (69-346AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Krueppel-like factor 5 (Basic transcription element-binding protein 2) (BTE-binding protein 2) (Colon krueppel-like factor) (GC-box-binding protein 2) (Intestinal-enriched krueppel-like factor) (Transcription factor BTEB2), KLF5, BTEB2 CKLF IKLF
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	KLF5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.