



KLF5 Antibody, FITC conjugated

Product Code	CSB-PA613429LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13887
Immunogen	Recombinant Human Krueppel-like factor 5 protein (69-346AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Krueppel-like factor 5 (Basic transcription element-binding protein 2) (BTE-binding protein 2) (Colon krueppel-like factor) (GC-box-binding protein 2) (Intestinal-enriched krueppel-like factor) (Transcription factor BTEB2), KLF5, BTEB2 CKLF IKLF
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	KLF5