



MLC1 Antibody, FITC conjugated

Product Code	CSB-PA613581LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15049
Immunogen	Recombinant Human Membrane protein MLC1 protein (1-49AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Membrane protein MLC1, MLC1, KIAA0027 WKL1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	MLC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.