



# MLC1 Antibody, Biotin conjugated

Product Code	CSB-PA613581LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15049
Immunogen	Recombinant Human Membrane protein MLC1 protein (1-49AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Membrane protein MLC1, MLC1, KIAA0027 WKL1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	MLC1