





ACOX1 Antibody

Product Code	CSB-PA617998HA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15067
Immunogen	Recombinant Human Peroxisomal acyl-coenzyme A oxidase 1 protein (113-273AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peroxisomal acyl-coenzyme A oxidase 1 (AOX) (EC 1.3.3.6) (Palmitoyl-CoA oxidase) (Straight-chain acyl-CoA oxidase) (SCOX), ACOX1, ACOX
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ACOX1


Image



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA617998HA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA617998HA01HU at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA617998HA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.