



PRKAA1 Antibody, Biotin conjugated

Product Code	CSB-PA618759LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13131
Immunogen	Recombinant Human 5'-AMP-activated protein kinase catalytic subunit alpha-1 protein (325-543AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	5'-AMP-activated protein kinase catalytic subunit alpha-1 (AMPK subunit alpha-1) (EC 2.7.11.1) (Acetyl-CoA carboxylase kinase) (ACACA kinase) (EC 2.7.11.27) (Hydroxymethylglutaryl-CoA reductase kinase) (HMGCR kinase) (EC 2.7.11.31) (Tau-protein kinase PRKAA1) (EC 2.7.11.26), PRKAA1, AMPK1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PRKAA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.