




EXT1 Antibody

Product Code	CSB-PA620959LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16394
Immunogen	Recombinant Human Exostosin-1 protein (28-171AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Exostosin-1 (EC 2.4.1.224) (EC 2.4.1.225) (Glucuronosyl-N-acetylglucosaminyl-proteoglycan/N-acetylglucosaminyl-proteoglycan 4-alpha-N-acetylglucosaminyltransferase) (Multiple exostoses protein 1) (Putative tumor suppressor protein EXT1), EXT1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	EXT1

Image



Western Blot

Positive WB detected in: HEK293 whole cell lysate

All lanes: EXT1 antibody at 3.3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 87 kDa

Observed band size: 87 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.