



EXT1 Antibody, FITC conjugated

Product Code	CSB-PA620959LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16394
Immunogen	Recombinant Human Exostosin-1 protein (28-171AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Exostosin-1 (EC 2.4.1.224) (EC 2.4.1.225) (Glucuronosyl-N-acetylglucosaminyl-proteoglycan/N-acetylglucosaminyl-proteoglycan 4-alpha-N-acetylglucosaminyltransferase) (Multiple exostoses protein 1) (Putative tumor suppressor protein EXT1), EXT1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	EXT1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.