



SWI5 Antibody, Biotin conjugated

Product Code	CSB-PA636328LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q1ZZU3
Immunogen	Recombinant Human DNA repair protein SWI5 homolog protein (54-112AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA repair protein SWI5 homolog (HBV DNAPTP1-transactivated protein A) (Protein SAE3 homolog), SWI5, C9orf119 SAE3
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	SWI5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.