




ZRANB3 Antibody

Product Code	CSB-PA705058LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5FWF4
Immunogen	Recombinant Human DNA annealing helicase and endonuclease ZRANB3 protein (370-624AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA annealing helicase and endonuclease ZRANB3 (Annealing helicase 2) (AH2) (Zinc finger Ran-binding domain-containing protein 3) [Includes: DNA annealing helicase ZRANB3 (EC 3.6.4.-); Endonuclease ZRANB3 (EC 3.1.-.-)], ZRANB3
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZRANB3

Image



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA705058LA01HU at dilution of 1:100

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.