



PREX2 Antibody, HRP conjugated

Product Code	CSB-PA741122LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q70Z35
Immunogen	Recombinant Human Phosphatidylinositol 3,4,5-trisphosphate-dependent Rac exchanger 2 protein (1417-1592AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphatidylinositol 3,4,5-trisphosphate-dependent Rac exchanger 2 protein (P-Rex2) (PtdIns(3,4,5)-dependent Rac exchanger 2) (DEP domain-containing protein 2), PREX2, DEPDC2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PREX2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.