



PLEKHA7 Antibody, FITC conjugated

Product Code	CSB-PA750339LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6IQ23
Immunogen	Recombinant Human Pleckstrin homology domain-containing family A member 7 protein (294-654AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pleckstrin homology domain-containing family A member 7 (PH domain-containing family A member 7), PLEKHA7
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PLEKHA7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.