





GPR179 Antibody

Product Code	CSB-PA757664LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6PRD1
Immunogen	Recombinant Human Probable G-protein coupled receptor 179 protein (1440-1671AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable G-protein coupled receptor 179 (Probable G-protein coupled receptor 158-like 1) (GPR158-like), GPR179, GPR158L GPR158L1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GPR179

Image



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA757664LA01HU at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA757664LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.