



GPR179 Antibody, FITC conjugated

Product Code	CSB-PA757664LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6PRD1
Immunogen	Recombinant Human Probable G-protein coupled receptor 179 protein (1440-1671AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable G-protein coupled receptor 179 (Probable G-protein coupled receptor 158-like 1) (GPR158-like), GPR179, GPR158L GPR158L1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GPR179
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.