



LINC00299 Antibody, HRP conjugated

Product Code	CSB-PA765077LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6ZSB3
Immunogen	Recombinant Human Putative uncharacterized protein encoded by LINC00299 protein (2-105AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative uncharacterized protein encoded by LINC00299, LINC00299, C2orf46 NCRNA00299
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	LINC00299
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.