



LINC00299 Antibody, Biotin conjugated

Product Code	CSB-PA765077LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6ZSB3
Immunogen	Recombinant Human Putative uncharacterized protein encoded by LINC00299 protein (2-105AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative uncharacterized protein encoded by LINC00299, LINC00299, C2orf46 NCRNA00299
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	LINC00299