



SLC51A Antibody, HRP conjugated

Product Code	CSB-PA768212LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86UW1
Immunogen	Recombinant Human Organic solute transporter subunit alpha protein (1-48AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Organic solute transporter subunit alpha (OST-alpha) (Solute carrier family 51 subunit alpha), SLC51A, OSTA
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC51A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.