



SLC51A Antibody, HRP conjugated

Product Code	CSB-PA768212LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86UW1
Immunogen	Recombinant Human Organic solute transporter subunit alpha protein (1-48AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Organic solute transporter subunit alpha (OST-alpha) (Solute carrier family 51 subunit alpha), SLC51A, OSTA
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC51A