



# SLC51B Antibody, FITC conjugated

Product Code	CSB-PA771440LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86UW2
Immunogen	Recombinant Human Organic solute transporter subunit beta protein (57-128AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Organic solute transporter subunit beta (OST-beta) (Solute carrier family 51 subunit beta), SLC51B, OSTB
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC51B