



SYVN1 Antibody, HRP conjugated

Product Code	CSB-PA803115LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86TM6
Immunogen	Recombinant Human E3 ubiquitin-protein ligase synoviolin protein (338-487AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase synoviolin (EC 2.3.2.27) (RING-type E3 ubiquitin transferase synoviolin) (Synovial apoptosis inhibitor 1), SYVN1, HRD1 KIAA1810
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SYVN1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.