



SIRT2 Antibody, HRP conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA812889LB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q8IXJ6 |
| Immunogen | Recombinant Human NAD-dependent protein deacetylase sirtuin-2 protein (177-279AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | NAD-dependent protein deacetylase sirtuin-2 (EC 3.5.1.-) (Regulatory protein SIR2 homolog 2) (SIR2-like protein 2), SIRT2, SIR2L SIR2L2 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Epigenetics and Nuclear Signaling |
| Target Names | SIRT2 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |