



SIRT2 Antibody, FITC conjugated

Product Code	CSB-PA812889LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IXJ6
Immunogen	Recombinant Human NAD-dependent protein deacetylase sirtuin-2 protein (177-279AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NAD-dependent protein deacetylase sirtuin-2 (EC 3.5.1.-) (Regulatory protein SIR2 homolog 2) (SIR2-like protein 2), SIRT2, SIR2L SIR2L2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	SIRT2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.