



WFDC8 Antibody, HRP conjugated

Product Code	CSB-PA818226LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IUA0
Immunogen	Recombinant Human WAP four-disulfide core domain protein 8 protein (44-94AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	WAP four-disulfide core domain protein 8 (Putative protease inhibitor WAP8), WFDC8, C20orf170 WAP8
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	WFDC8