



WFDC8 Antibody, FITC conjugated

Product Code	CSB-PA818226LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IUA0
Immunogen	Recombinant Human WAP four-disulfide core domain protein 8 protein (44-94AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	WAP four-disulfide core domain protein 8 (Putative protease inhibitor WAP8), WFDC8, C20orf170 WAP8
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	WFDC8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.