





SLC30A8 Antibody

Product Code	CSB-PA818247LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IWU4
Immunogen	Recombinant Human Zinc transporter 8 protein (320-369AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc transporter 8 (ZnT-8) (Solute carrier family 30 member 8), SLC30A8, ZNT8
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC30A8

Image



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA818247LA01HU at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA818247LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)