



SLC30A8 Antibody, FITC conjugated

Product Code	CSB-PA818247LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IWU4
Immunogen	Recombinant Human Zinc transporter 8 protein (320-369AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc transporter 8 (ZnT-8) (Solute carrier family 30 member 8), SLC30A8, ZNT8
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC30A8