



# DOCK7 Antibody, HRP conjugated

|                     |  |
|---------------------|--|
| Product Code        | CSB-PA836266LB01HU   |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                    |
| Uniprot No.         | Q96N67   |
| Immunogen           | Recombinant Human Dedicator of cytokinesis protein 7 protein (1440-1671AA)       |
| Raised In           | Rabbit   |
| Species Reactivity  | Human  |
| Tested Applications | ELISA  |
| Form                | Liquid   |
| Conjugate           | HRP  |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified   |
| Isotype             | IgG  |
| Clonality           | Polyclonal   |
| Alias               | Dedicator of cytokinesis protein 7, DOCK7, KIAA1771                              |
| Immunogen Species   | Homo sapiens (Human)   |
| Research Area       | Signal Transduction  |
| Target Names        | DOCK7  |