




SLC30A7 Antibody

Product Code	CSB-PA836732LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NEW0
Immunogen	Recombinant Human Zinc transporter 7 protein (162-236AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc transporter 7 (ZnT-7) (Solute carrier family 30 member 7) (Znt-like transporter 2), SLC30A7, ZNT7 ZNTL2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC30A7

Image



Western Blot

Positive WB detected in: Mouse heart tissue,
Mouse skeletal muscle tissue


All lanes: SLC30A7 antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 42 kDa

Observed band size: 42 kDa



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA836732LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.